

NOS3 / eNOS Rabbit anti-Human Polyclonal (Ser1177) (PE) Antibody - LS-C277121 - LSBio	
CatalogID:	LS-C277121
Target:	nitric oxide synthase 3 (endothelial cell) (NOS3)
Synonyms:	NOS3 Antibody, Constitutive NOS Antibody, CNOS Antibody, Endothelial NOS Antibody, ECNOS Antibody, NOS type III Antibody, NOSIII Antibody, EC-NOS Antibody, ENOS Antibody, Nitric oxide synthase 3 Antibody
Host	NOS3 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	NOS3 / eNOS antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	NOS3 / eNOS antibody was raised against kLH-conjugated synthetic phosphopeptide from amino acids around Ser1177 of human eNos.
Specificity:	Human NOS3 / eNOS
Epitope:	Ser1177
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/30/2014	
© 2014 LifeSpan BioSciences	