

ERBB4 / HER4 Rabbit anti-Human Polyclonal (Tyr1188) (PE) Antibody - LS-C277168 - LSBio	
CatalogID:	LS-C277168
Target:	v-erb-b2 avian erythroblastic leukemia viral oncogene homolog 4 (ERBB4)
Synonyms:	ERBB4 Antibody, ERBB 4 Antibody, HER4 Antibody, p180 erb4 Antibody, p180erbB4 Antibody
Family / Subfamily:	Protein Kinase / EGF Receptor
Host	ERBB4 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	ERBB4 / HER4 antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	ERBB4 / HER4 antibody was raised against kLH-conjugated synthetic phosphopeptide from amino acids around Tyr1188 of human HER4.
Specificity:	Human ERBB4 / HER4
Epitope:	Tyr1188
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only Not for resale without prior written consent from LifeSpan BioSciences, Inc. Created on 9/30/2014 © 2014 LifeSpan BioSciences	