

| ESR2 / ER Beta Rabbit anti-Human Polyclonal (aa376-405) (FITC) Antibody - LS-C277183 - LSBio |  |
|--|--|
| <b>CatalogID:</b>  | LS-C277183   |
| <b>Target:</b>   | estrogen receptor 2 (ER beta) (ESR2)   |
| <b>Synonyms:</b>   | ESR2 Antibody, ESTRB Antibody, ER Beta Antibody, ERbeta Antibody, ER-BETA Antibody, ESR-BETA Antibody, ESRB Antibody, Estrogen receptor 2 (ER beta) Antibody, NR3A2 Antibody, Estrogen receptor beta Antibody, Estrogen receptor beta 4 Antibody |
| <b>Family / Subfamily:</b>   | NHR / NR3 Steroid receptor   |
| <b>Host</b>  | ESR2 antibody was produced in Rabbit   |
| <b>Clonality:</b>  | Polyclonal   |
| <b>Isotype:</b>  | IgG  |
| <b>Conjugations:</b>   | Fluorescein (FITC)   |
| <b>Immunogen Species:</b>  | ESR2 / ER Beta antibody was raised against Human   |
| <b>Antigen Type:</b>   | Synthetic peptide  |
| <b>Immunogen:</b>  | ESR2 / ER Beta antibody was raised against kLH-conjugated synthetic peptide (aa376-405) from internal region of human ESR2.  |
| <b>Specificity:</b>  | Human ER Beta/Estrogen Receptor Beta   |
| <b>Epitope:</b>  | aa376-405  |
| <b>Reactivity:</b>   | Human, Mouse   |
| <b>Purification:</b>   | Affinity purified  |
| <b>Presentation:</b>   | PBS, pH 7.2, 0.09% sodium azide.   |
| <b>Recommended Storage:</b>  | May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.   |
| <b>Usage Summary:</b>  | The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.   |
| <b>Uses:</b>   | IHC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)   |
| <b>Size:</b>   | 200 µl   |
| <b>Requested From:</b>   | Japan  |
| <b>Laboratory Reagent For In Vitro Research Use Only</b>                                     |  |
| <b>Not for resale without prior written consent from LifeSpan BioSciences, Inc.</b>          |  |
| Created on 9/30/2014<br>© 2014 LifeSpan BioSciences  |  |