

FGR Rabbit anti-Human Polyclonal (N-Terminus) (PE) Antibody - LS-C277340 - LSBio

CatalogID:	LS-C277340
Target:	feline Gardner-Rasheed sarcoma viral oncogene homolog (FGR)
Synonyms:	FGR Antibody, C-fgr protooncogene Antibody, C-src2 Antibody, C-fgr Antibody, Fgr kinase Antibody, Fgr oncogene Antibody, p55-Fgr Antibody, p58-Fgr Antibody, p58c-fgr Antibody, Proto-oncogene c-Fgr Antibody, Tyrosine-protein kinase Fgr Antibody, p55c-fgr Antibody, SRC2 Antibody, C-src-2 proto-oncogene Antibody, Fgr proto oncogene Antibody, p55-c-fgr protein Antibody
Family / Subfamily:	Protein Kinase / Src
Host	FGR antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	FGR antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	FGR antibody was raised against kLH-conjugated synthetic peptide from the N-terminal region of human FGR.
Specificity:	Human FGR.
Epitope:	N-Terminus
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences