

GAK Rabbit anti-Human Polyclonal (N-Terminus) (PE) Antibody - LS-C277489 - LSBio

CatalogID:	LS-C277489
Target:	cyclin G associated kinase (GAK)
Synonyms:	GAK Antibody, Auxilin-2 Antibody, Cyclin-G-associated kinase Antibody, Cyclin G associated kinase Antibody, DNAJ26 Antibody, DNAJC26 Antibody
Family / Subfamily:	Protein Kinase / NAK
Host	GAK antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	GAK antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	GAK antibody was raised against kLH-conjugated synthetic peptide from the N-terminal region of human GAK.
Specificity:	Human GAK
Epitope:	N-Terminus
Reactivity:	Human, Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences