

PIGK Rabbit anti-Human Polyclonal (N-Terminus) (FITC) Antibody - LS-C277653 - LSBio	
CatalogID:	LS-C277653
Target:	phosphatidylinositol glycan anchor biosynthesis, class K (PIGK)
Synonyms:	PIGK Antibody, GPI transamidase Antibody, GPI transamidase subunit Antibody, GPI-anchor transamidase Antibody, GPI8 Antibody, GPI8 homolog Antibody, HGPI8 Antibody, PIG-K Antibody
Family / Subfamily:	Protease / Cysteine C13
Host	PIGK antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	PIGK antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	PIGK antibody was raised against kLH-conjugated synthetic peptide from the N-terminal region of human GPI8.
Specificity:	Human PIGK
Epitope:	N-Terminus
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/30/2014	
© 2014 LifeSpan BioSciences	