

HGF Rabbit anti-Human Polyclonal (aa529-559) (FITC) Antibody - LS-C277994 - LSBio

CatalogID:	LS-C277994
Target:	hepatocyte growth factor (hepapoietin A; scatter factor) (HGF)
Synonyms:	HGF Antibody, DFNB39 Antibody, Hepatocyte growth factor Antibody, Hepatopoeitin-A Antibody, Hepapoietin A Antibody, Hepatopoeitin-A Antibody, Scatter factor) Antibody, SF Antibody, F-TCF Antibody, HGFB Antibody, HPTA Antibody, Scatter factor Antibody
Family / Subfamily:	Protease / Serine S1
Host	HGF antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	HGF antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	HGF antibody was raised against kLH-conjugated synthetic peptide (aa529-559) from the C-terminus of human HGF.
Specificity:	Human HGF
Epitope:	aa529-559
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences