

HK1 / Hexokinase 1 Rabbit anti-Human Polyclonal (C-Terminus) (FITC) Antibody - LS-C278030 - LSBio

CatalogID:	LS-C278030
Target:	hexokinase 1 (HK1)
Synonyms:	HK1 Antibody, Brain form hexokinase Antibody, Hexokinase 1 Antibody, Hexokinase type I Antibody, HK I Antibody, HK1-tb Antibody, HK1-tc Antibody, HKI Antibody, Hexokinase type 1 Antibody, HK1-ta Antibody, HXK1 Antibody, Glycolytic enzyme Antibody, Hexokinase-1 Antibody, Type 1 hexokinase Antibody, Type i hexokinase Antibody
Family / Subfamily:	Non-protein Kinase / Hexokinase
Host	HK1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	HK1 / Hexokinase 1 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	HK1 / Hexokinase 1 antibody was raised against kLH-conjugated synthetic peptide from the C-terminus of human HK1.
Specificity:	Human HK1 / Hexokinase 1
Epitope:	C-Terminus
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences