

**HEXA Rabbit anti-Human Polyclonal (Internal) (FITC) Antibody - LS-C278084 - LSBio**

|                             |  |
|-----------------------------|--|
| <b>CatalogID:</b>           | LS-C278084   |
| <b>Target:</b>              | hexosaminidase A (alpha polypeptide) (HEXA)  |
| <b>Synonyms:</b>            | HEXA Antibody, Hexosaminidase subunit A Antibody, TSD Antibody   |
| <b>Host</b>                 | HEXA antibody was produced in Rabbit   |
| <b>Clonality:</b>           | Polyclonal   |
| <b>Isotype:</b>             | IgG  |
| <b>Conjugations:</b>        | Fluorescein (FITC)   |
| <b>Immunogen Species:</b>   | HEXA antibody was raised against Human   |
| <b>Antigen Type:</b>        | Synthetic peptide  |
| <b>Immunogen:</b>           | HEXA antibody was raised against kLH-conjugated synthetic peptide from internal region of human HEXA.  |
| <b>Specificity:</b>         | Human HEXA.  |
| <b>Epitope:</b>             | Internal   |
| <b>Reactivity:</b>          | Human  |
| <b>Purification:</b>        | Affinity purified  |
| <b>Presentation:</b>        | PBS, pH 7.2, 0.09% sodium azide.   |
| <b>Recommended Storage:</b> | May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months. |
| <b>Usage Summary:</b>       | The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.                           |
| <b>Uses:</b>                | Western blot, ELISA (Optimal dilution to be determined by the researcher)  |
| <b>Size:</b>                | 200 µl   |
| <b>Requested From:</b>      | Japan  |

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences