

ITPKA Rabbit anti-Human Polyclonal (C-Terminus) (FITC) Antibody - LS-C278334 - LSBio	
CatalogID:	LS-C278334
Target:	inositol-trisphosphate 3-kinase A (ITPKA)
Synonyms:	ITPKA Antibody, IP3-3KA Antibody, IP3KA Antibody, InsP 3-kinase A Antibody, IP3 3-kinase A Antibody, IP3K A Antibody
Family / Subfamily:	Non-protein Kinase / ITPK
Host	ITPKA antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	ITPKA antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	ITPKA antibody was raised against kLH-conjugated synthetic peptide from the C-terminus of human ITPKA.
Specificity:	Human ITPKA
Epitope:	C-Terminus
Reactivity:	Human, Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only Not for resale without prior written consent from LifeSpan BioSciences, Inc. Created on 9/30/2014 © 2014 LifeSpan BioSciences	