

IL1A / IL1 Alpha Rabbit anti-Human Polyclonal (Ser87) (PE) Antibody - LS-C278471 - LSBio	
<b>CatalogID:</b>	LS-C278471
<b>Target:</b>	interleukin 1, alpha (IL1A)
<b>Synonyms:</b>	IL1A Antibody, Hematopoietin-1 Antibody, IL-1 alpha Antibody, IL1 Antibody, IL1-ALPHA Antibody, IL1 Alpha Antibody, Interleukin-1 alpha Antibody, IL-1A Antibody, IL1F1 Antibody, Pro-interleukin-1-alpha Antibody, Preinterleukin 1 alpha Antibody, Interleukin 1, alpha Antibody
<b>Family / Subfamily:</b>	IL-1 / not assigned-IL-1
<b>Host</b>	IL1A antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Phycoerythrin (PE)
<b>Immunogen Species:</b>	IL1A / IL1 Alpha antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	IL1A / IL1 Alpha antibody was raised against kLH-conjugated synthetic peptide from the amino acid residues around Ser87 of human IL1A.
<b>Specificity:</b>	Human IL1A / IL1 Alpha.
<b>Epitope:</b>	Ser87
<b>Reactivity:</b>	Human
<b>Purification:</b>	Affinity purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	Western blot, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/30/2014	
© 2014 LifeSpan BioSciences	