

IL1A / IL1 Alpha Rabbit anti-Human Polyclonal (Internal) (PE) Antibody - LS-C278475 - LSBio	
CatalogID:	LS-C278475
Target:	interleukin 1, alpha (IL1A)
Synonyms:	IL1A Antibody, Hematopoietin-1 Antibody, IL-1 alpha Antibody, IL1 Antibody, IL1-ALPHA Antibody, IL1 Alpha Antibody, Interleukin-1 alpha Antibody, IL-1A Antibody, IL1F1 Antibody, Pro-interleukin-1-alpha Antibody, Preinterleukin 1 alpha Antibody, Interleukin 1, alpha Antibody
Family / Subfamily:	IL-1 / not assigned-IL-1
Host	IL1A antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	IL1A / IL1 Alpha antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	IL1A / IL1 Alpha antibody was raised against KLH-conjugated synthetic peptide from internal region of human IL1A.
Specificity:	Human IL1A / IL1 Alpha.
Epitope:	Internal
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/30/2014	
© 2014 LifeSpan BioSciences	