

JUN / c-Jun Rabbit anti-Human Polyclonal (Ser63) (PE) Antibody - LS-C278649 - LSBio	
CatalogID:	LS-C278649
Target:	jun proto-oncogene
Synonyms:	JUN Antibody, Activator protein 1 Antibody, AP-1 Antibody, C-Jun Antibody, Jun proto-oncogene Antibody, p39 Antibody, Jun oncogene Antibody, Proto-oncogene c-Jun Antibody, AP1 Antibody, Enhancer-binding protein AP1 Antibody, G0s7 Antibody, Transcription factor AP-1 Antibody, v-jun Antibody
Family / Subfamily:	Transcription factor / not assigned-Transcription factor
Host	JUN antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	JUN / c-Jun antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	JUN / c-Jun antibody was raised against kLH-conjugated synthetic phosphopeptide from amino acids around Ser63 of human cJun.
Specificity:	Human JUN / c-Jun
Epitope:	Ser63
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/30/2014	
© 2014 LifeSpan BioSciences	