

LRP1 / CD91 Rabbit anti-Human Polyclonal (Ser452) (PE) Antibody - LS-C278947 - LSBio

CatalogID:	LS-C278947
Target:	low density lipoprotein receptor-related protein 1 (LRP1)
Synonyms:	LRP1 Antibody, Apolipoprotein E receptor Antibody, CD91 antigen Antibody, CD91 Antibody, A2MR Antibody, IGFBP3R Antibody, Alpha-2-macroglobulin receptor Antibody, TGFBR5 Antibody, Type V tgf-beta receptor Antibody, APOER Antibody, APR Antibody, LRP Antibody, LRP-1 Antibody
Host	LRP1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	LRP1 / CD91 antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	LRP1 / CD91 antibody was raised against kLH-conjugated synthetic phosphopeptide from amino acids around Ser452 of human LPR1.
Specificity:	Human LRP1 / LRP
Epitope:	Ser452
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences