

MAPKAPK2 / MK2 Rabbit anti-Human Polyclonal (Ser272) (PE) Antibody - LS-C279210 - LSBio	
CatalogID:	LS-C279210
Target:	mitogen-activated protein kinase-activated protein kinase 2 (MAPKAPK2)
Synonyms:	MAPKAPK2 Antibody, MK-2 Antibody, MAP KAP Kinase 2 Antibody, MAPKAP kinase 2 Antibody, MK2 Antibody, MAPKAP-K2 Antibody, MAPKAPK-2 Antibody
Family / Subfamily:	Protein Kinase / MAPKAPK
Host	MAPKAPK2 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	MAPKAPK2 / MK2 antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	MAPKAPK2 / MK2 antibody was raised against kLH-conjugated synthetic phosphopeptide from amino acids around Ser272 of human MAPKAPK2.
Specificity:	Human MAPKAPK2 / MK2
Epitope:	Ser272
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/30/2014	
© 2014 LifeSpan BioSciences	