

MMP3 / MMP-3 Rabbit anti-Human Polyclonal (Internal) (FITC) Antibody - LS-C279278 - LSBio	
CatalogID:	LS-C279278
Target:	matrix metalloproteinase 3 (stromelysin 1, progelatinase) (MMP3)
Synonyms:	MMP3 Antibody, CHDS6 Antibody, MMP-3 Antibody, Prostromelysin-1 Antibody, Proteoglycanase Antibody, SL-1 Antibody, STR1 Antibody, Matrix metalloproteinase-3 Antibody, Transin Antibody, Transin-1 Antibody, SLN-1 Antibody, STMY Antibody, Progelatinase Antibody, STMY1 Antibody, Stromelysin-1 Antibody
Family / Subfamily:	Protease / Metalloproteinase M10A
Host	MMP3 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	MMP3 / MMP-3 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	MMP3 / MMP-3 antibody was raised against kLH-conjugated synthetic peptide from internal region of human MMP3.
Specificity:	Human MMP3 / MMP-3
Epitope:	Internal
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/30/2014 © 2014 LifeSpan BioSciences	