

MMP16 / MT3-MMP Rabbit anti-Human Polyclonal (N-Terminus) (PE) Antibody - LS-C279319 - LSBio

CatalogID:	LS-C279319
Target:	matrix metalloproteinase 16 (membrane-inserted) (MMP16)
Synonyms:	MMP16 Antibody, C8orf57 Antibody, MMP-X2 Antibody, MT-MMP3 Antibody, MT3MMP Antibody, MT-MMP 3 Antibody, MT-MMP2 Antibody, MT2-MMP Antibody, Matrix metalloproteinase-16 Antibody, Membrane type 3-MMP Antibody, MMP-16 Antibody, MMPX2 Antibody, MT3-MMP Antibody, MTMMP3 Antibody
Family / Subfamily:	Protease / Metalloproteinase M10A
Host	MMP16 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	MMP16 / MT3-MMP antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	MMP16 / MT3-MMP antibody was raised against kLH-conjugated synthetic peptide from the N-terminal region of human MMP16.
Specificity:	Human MMP16 / MT3-MMP
Epitope:	N-Terminus
Reactivity:	Human, Hamster
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences