

MELK Rabbit anti-Human Polyclonal (C-Terminus) (FITC) Antibody - LS-C279443 - LSBio	
CatalogID:	LS-C279443
Target:	maternal embryonic leucine zipper kinase (MELK)
Synonyms:	MELK Antibody, KIAA0175 Antibody, HPK38 Antibody, Protein kinase Eg3 Antibody, Tyrosine-protein kinase MELK Antibody, Protein kinase PK38 Antibody, HMELK Antibody, Mpk38 Antibody, PEg3 kinase Antibody
Family / Subfamily:	Protein Kinase / CAMKL
Host	MELK antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	MELK antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	MELK antibody was raised against kLH-conjugated synthetic peptide from the C-terminus of human MELK.
Specificity:	Human MELK.
Epitope:	C-Terminus
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/30/2014	
© 2014 LifeSpan BioSciences	