

MOS Rabbit anti-Human Polyclonal (N-Terminus) (PE) Antibody - LS-C279656 - LSBio

CatalogID:	LS-C279656
Target:	v-mos Moloney murine sarcoma viral oncogene homolog
Synonyms:	MOS Antibody, Oocyte maturation factor mos Antibody, Proto-oncogene c-Mos Antibody, MSV Antibody, C-mos Antibody
Family / Subfamily:	Protein Kinase / MOS
Host	MOS antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	MOS antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	MOS antibody was raised against kLH-conjugated synthetic peptide from the N-terminal region of human MOS.
Specificity:	Human MOS
Epitope:	N-Terminus
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences