

MEF2C Rabbit anti-Human Polyclonal (Ser396) (PE) Antibody - LS-C279811 - LSBio

CatalogID:	LS-C279811
Target:	myocyte enhancer factor 2C (MEF2C)
Synonyms:	MEF2C Antibody, C5DELq14.3 Antibody, DEL5q14.3 Antibody, Myocyte enhancer factor 2C Antibody
Host	MEF2C antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	MEF2C antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	MEF2C antibody was raised against kLH-conjugated synthetic phosphopeptide from amino acids around Ser396 of human MEF2C.
Specificity:	Human MEF2C.
Epitope:	Ser396
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences