

PKMYT1 / MYT1 Rabbit anti-Human Polyclonal (Thr495) (FITC) Antibody - LS-C279854 - LSBio	
CatalogID:	LS-C279854
Target:	protein kinase, membrane associated tyrosine/threonine 1 (PKMYT1)
Synonyms:	PKMYT1 Antibody, Myt1 kinase Antibody
Family / Subfamily:	Protein Kinase / wee1/mik1
Host	PKMYT1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	PKMYT1 / MYT1 antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	PKMYT1 / MYT1 antibody was raised against kLH-conjugated synthetic phosphopeptide from amino acids around Thr495 of human MYT1.
Specificity:	Human PKMYT1 / MYT1.
Epitope:	Thr495
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/30/2014	
© 2014 LifeSpan BioSciences	