

NANOG Rabbit anti-Human Polyclonal (Ser71) (FITC) Antibody - LS-C279945 - LSBio

CatalogID:	LS-C279945
Target:	Nanog homeobox
Synonyms:	NANOG Antibody, Nanog homeobox Antibody, NANOGP8 Antibody, HNanog Antibody, Homeobox protein NANOG Antibody
Host	NANOG antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	NANOG antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	NANOG antibody was raised against kLH-conjugated synthetic phosphopeptide from amino acids around Ser71 of human NANOG.
Specificity:	Human NANOG.
Epitope:	Ser71
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences