

NANOG Rabbit anti-Human Polyclonal (Ser285) (Azide-free) (HRP) Antibody - LS-C279951 - LSBio	
CatalogID:	LS-C279951
Target:	Nanog homeobox
Synonyms:	NANOG Antibody, Nanog homeobox Antibody, NANOGP8 Antibody, HNanog Antibody, Homeobox protein NANOG Antibody
Host	NANOG antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Horseradish Peroxidase (HRP)
Modifications:	Azide-free
Immunogen Species:	NANOG antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	NANOG antibody was raised against kLH-conjugated synthetic phosphopeptide from amino acids around Ser285 of human NANOG.
Specificity:	Human NANOG.
Epitope:	Ser285
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, azide-free.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/30/2014	
© 2014 LifeSpan BioSciences	