

NDRG1 Rabbit anti-Human Polyclonal (N-Terminus) (Azide-free) (HRP) Antibody - LS-C280155 - LSBio

CatalogID:	LS-C280155
Target:	N-myc downstream regulated 1 (NDRG1)
Synonyms:	NDRG1 Antibody, CMT4D Antibody, DRG-1 Antibody, GC4 Antibody, HMSNL Antibody, Protein regulated by oxygen-1 Antibody, PROXY1 Antibody, Protein NDRG1 Antibody, TARG1 Antibody, TDD5 Antibody, NDR1 Antibody, RTP Antibody, CAP43 Antibody, N-myc downstream regulated 1 Antibody, NMSL Antibody, RIT42 Antibody
Host	NDRG1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Horseradish Peroxidase (HRP)
Modifications:	Azide-free
Immunogen Species:	NDRG1 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	NDRG1 antibody was raised against kLH-conjugated synthetic peptide from the N-terminal region of human NDRG1.
Specificity:	Human NDRG1
Epitope:	N-Terminus
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, azide-free.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences