

**NDRG1 Rabbit anti-Human Polyclonal (N-Terminus) (PE) Antibody - LS-C280156 - LSBio**

<b>CatalogID:</b>	LS-C280156
<b>Target:</b>	N-myc downstream regulated 1 (NDRG1)
<b>Synonyms:</b>	NDRG1 Antibody, CMT4D Antibody, DRG-1 Antibody, GC4 Antibody, HMSNL Antibody, Protein regulated by oxygen-1 Antibody, PROXY1 Antibody, Protein NDRG1 Antibody, TARG1 Antibody, TDD5 Antibody, NDR1 Antibody, RTP Antibody, CAP43 Antibody, N-myc downstream regulated 1 Antibody, NMSL Antibody, RIT42 Antibody
<b>Host</b>	NDRG1 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Phycoerythrin (PE)
<b>Immunogen Species:</b>	NDRG1 antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	NDRG1 antibody was raised against kLH-conjugated synthetic peptide from the N-terminal region of human NDRG1.
<b>Specificity:</b>	Human NDRG1
<b>Epitope:</b>	N-Terminus
<b>Reactivity:</b>	Human
<b>Purification:</b>	Affinity purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences