

NR0B1 / DAX1 Rabbit anti-Human Polyclonal (N-Terminus) (FITC) Antibody - LS-C280268 - LSBio	
<b>CatalogID:</b>	LS-C280268
<b>Target:</b>	nuclear receptor subfamily 0, group B, member 1 (NR0B1)
<b>Synonyms:</b>	NR0B1 Antibody, AHX Antibody, AHC Antibody, AHCH Antibody, DAX1 Antibody, DSS Antibody, DAX 1 Antibody, DAX-1 Antibody, Adrenal hypoplasia congenita Antibody, Dosage-sensitive sex reversal Antibody, GTD Antibody, HHG Antibody, Nuclear hormone receptor Antibody, NROB1 Antibody, Nuclear receptor DAX-1 Antibody, SRXY2 Antibody, Nuclear receptor DAX1 Antibody, Dax-1 receptor Antibody
<b>Family / Subfamily:</b>	NHR / NR0 Knirps-like
<b>Host</b>	NR0B1 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	NR0B1 / DAX1 antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	NR0B1 / DAX1 antibody was raised against kLH-conjugated synthetic peptide from the N-terminal region of human NROB1.
<b>Specificity:</b>	Human NR0B1 / DAX1
<b>Epitope:</b>	N-Terminus
<b>Reactivity:</b>	Human
<b>Purification:</b>	Affinity purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	Western blot, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/30/2014	
© 2014 LifeSpan BioSciences	