

OASL Rabbit anti-Human Polyclonal (C-Terminus) (PE) Antibody - LS-C280370 - LSBio

CatalogID:	LS-C280370
Target:	2'-5'-oligoadenylate synthetase-like (OASL)
Synonyms:	OASL Antibody, 2-5-OAS-related protein Antibody, 2-5-OAS-RP Antibody, p59-OASL Antibody, p59OASL Antibody, TR-interacting protein 14 Antibody, TRIP-14 Antibody, TRIP14 Antibody, OASLd Antibody, p59 OASL Antibody
Host	OASL antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	OASL antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	OASL antibody was raised against kLH-conjugated synthetic peptide from the C-terminus of human OASL.
Specificity:	Human OASL
Epitope:	C-Terminus
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences