

OSR1 Rabbit anti-Human Polyclonal (Thr310) (PE) Antibody - LS-C280438 - LSBio

CatalogID:	LS-C280438
Target:	odd-skipped related transcription factor 1 (OSR1)
Synonyms:	OSR1 Antibody, ODD Antibody, Protein odd-skipped-related 1 Antibody, Odd-skipped homolog Antibody
Host	OSR1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	OSR1 antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	OSR1 antibody was raised against kLH-conjugated synthetic phosphopeptide from amino acids around Thr310 of human OSR1.
Specificity:	Human OSR1
Epitope:	Thr310
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences