

POSTN / Periostin Rabbit anti-Human Polyclonal (aa694-724) (PE) Antibody - LS-C280444 - LSBio	
CatalogID:	LS-C280444
Target:	periostin, osteoblast specific factor (POSTN)
Synonyms:	POSTN Antibody, OSF2 Antibody, PDLPOSTN Antibody, Periostin isoform thy6 Antibody, Periostin isoform thy8 Antibody, RP11-412K4.1 Antibody, OSF-2 Antibody, Osteoblast-specific factor 2 Antibody, Periostin Antibody, Periostin isoform thy2 Antibody, Periostin isoform thy4 Antibody, PN Antibody
Host	POSTN antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	POSTN / Periostin antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	POSTN / Periostin antibody was raised against kLH-conjugated synthetic peptide (aa694-724) from the C-terminus of human Periostin.
Specificity:	Human POSTN / Periostin
Epitope:	aa694-724
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ICC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/30/2014	
© 2014 LifeSpan BioSciences	