

PDX1 / PDX-1 Rabbit anti-Human Polyclonal (Thr11) (FITC) Antibody - LS-C280690 - LSBio

CatalogID:	LS-C280690
Target:	pancreatic and duodenal homeobox 1 (PDX1)
Synonyms:	PDX1 Antibody, Glucose-sensitive factor Antibody, Insulin upstream factor 1 Antibody, IPF1 Antibody, Insulin promoter factor 1 Antibody, IPF-1 Antibody, IUF-1 Antibody, STF-1 Antibody, GSF Antibody, IDX-1 Antibody, Islet/duodenum homeobox-1 Antibody, IUF1 Antibody, MODY4 Antibody, PDX-1 Antibody
Host	PDX1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	PDX1 / PDX-1 antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	PDX1 / PDX-1 antibody was raised against kLH-conjugated synthetic phosphopeptide from amino acids around Thr11 of human PDX1.
Specificity:	Human PDX1 / PDX-1
Epitope:	Thr11
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences