

**PDGFRA / PDGFR Alpha Rabbit anti-Human Polyclonal (C-Terminus) (FITC) Antibody - LS-C281041 - LSBio**

<b>CatalogID:</b>	LS-C281041
<b>Target:</b>	platelet-derived growth factor receptor, alpha polypeptide (PDGFRA)
<b>Synonyms:</b>	PDGFRA Antibody, CD140A Antibody, CD140a antigen Antibody, PDGFalphaR Antibody, PDGFR2 Antibody, PDGFRalpha Antibody, AlphaPDGFR Antibody, RHEPDGFRA Antibody, PDGF receptor alpha Antibody, PDGFR-2 Antibody, PDGFRA/BCR fusion Antibody, APDGFR Antibody, Pdgf alpha-receptor Antibody, PDGF-R-alpha Antibody, PDGFR-alpha Antibody
<b>Family / Subfamily:</b>	Protein Kinase / PDGF Receptor
<b>Host</b>	PDGFRA antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	PDGFRA / PDGFR Alpha antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	PDGFRA / PDGFR Alpha antibody was raised against kLH-conjugated synthetic peptide from the C-terminus of human PDGFRA.
<b>Specificity:</b>	Human PDGFRA / PDGFR Alpha
<b>Epitope:</b>	C-Terminus
<b>Reactivity:</b>	Human, Mouse
<b>Purification:</b>	Affinity purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences