

| PDGFRA / PDGFR Alpha Rabbit anti-Human Polyclonal (Tyr720) (FITC) Antibody - LS-C281048 - LSBio | |
|---|---|
| CatalogID: | LS-C281048 |
| Target: | platelet-derived growth factor receptor, alpha polypeptide (PDGFRA) |
| Synonyms: | PDGFRA Antibody, CD140A Antibody, CD140a antigen Antibody, PDGFalphaR Antibody, PDGFR2 Antibody, PDGFRalpha Antibody, AlphaPDGFR Antibody, RHEPDGFRA Antibody, PDGF receptor alpha Antibody, PDGFR-2 Antibody, PDGFRA/BCR fusion Antibody, APDGFR Antibody, Pdgf alpha-receptor Antibody, PDGF-R-alpha Antibody, PDGFR-alpha Antibody |
| Family / Subfamily: | Protein Kinase / PDGF Receptor |
| Host | PDGFRA antibody was produced in Rabbit |
| Clonality: | Polyclonal |
| Isotype: | IgG |
| Conjugations: | Fluorescein (FITC) |
| Immunogen Species: | PDGFRA / PDGFR Alpha antibody was raised against Human |
| Antigen Type: | Synthetic peptide |
| Immunogen: | PDGFRA / PDGFR Alpha antibody was raised against kLH-conjugated synthetic peptide around Tyr720 of human PDGFRA. |
| Specificity: | Human PDGFRA / PDGFR Alpha |
| Epitope: | Tyr720 |
| Reactivity: | Human |
| Purification: | Affinity purified |
| Presentation: | PBS, pH 7.2, 0.09% sodium azide. |
| Recommended Storage: | May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months. |
| Usage Summary: | The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested. |
| Uses: | IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher) |
| Size: | 200 µl |
| Requested From: | Japan |
| Laboratory Reagent For In Vitro Research Use Only | |
| Not for resale without prior written consent from LifeSpan BioSciences, Inc. | |
| Created on 9/30/2014 | |
| © 2014 LifeSpan BioSciences | |