

PKLR Rabbit anti-Human Polyclonal (C-Terminus) (FITC) Antibody - LS-C281757 - LSBio	
CatalogID:	LS-C281757
Target:	pyruvate kinase, liver and RBC (PKLR)
Synonyms:	PKLR Antibody, PKL Antibody, Pyruvate kinase, liver and RBC Antibody, Pyruvate kinase type L Antibody, Red cell/liver pyruvate kinase Antibody, RPK Antibody, PK1 Antibody, PKR Antibody, Pyruvate kinase 1 Antibody, Pyruvate kinase isozymes R/L Antibody, R-type/L-type pyruvate kinase Antibody, PKRL Antibody, Pyruvate kinase isozyme R/L Antibody
Host	PKLR antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	PKLR antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	PKLR antibody was raised against kLH-conjugated synthetic peptide from the C-terminus of human PKLR.
Specificity:	Human PKLR
Epitope:	C-Terminus
Reactivity:	Human, Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/30/2014	
© 2014 LifeSpan BioSciences	