

**QKI Rabbit anti-Human Polyclonal (N-Terminus) (FITC) Antibody - LS-C281781 - LSBio**

<b>CatalogID:</b>	LS-C281781
<b>Target:</b>	QKI, KH domain containing, RNA binding
<b>Synonyms:</b>	QKI Antibody, Hqkl Antibody, HKQ Antibody, Protein quaking Antibody, QK3 Antibody, RNA binding protein HQK Antibody, Hqk Antibody, QK Antibody, QK1 Antibody
<b>Host</b>	QKI antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	QKI antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	QKI antibody was raised against kLH-conjugated synthetic peptide from the N-terminal region of human QKI.
<b>Specificity:</b>	Human QKI
<b>Epitope:</b>	N-Terminus
<b>Reactivity:</b>	Human
<b>Purification:</b>	Affinity purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences