

RAD9A / RAD9 Rabbit anti-Human Polyclonal (Ser328) (PE) Antibody - LS-C281896 - LSBio

CatalogID:	LS-C281896
Target:	RAD9 homolog A (S. pombe) (RAD9A)
Synonyms:	RAD9A Antibody, HRAD9 Antibody, RAD9 Antibody, RAD9 homolog Antibody, RAD9 (S. pombe) homolog Antibody, RAD9 homolog A (S. pombe) Antibody
Host	RAD9A antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	RAD9A / RAD9 antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	RAD9A / RAD9 antibody was raised against kLH-conjugated synthetic phosphopeptide from amino acids around Ser328 of human Rad9.
Specificity:	Human RAD9A / RAD9.
Epitope:	Ser328
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences