

RAF1 / RAF Rabbit anti-Human Polyclonal (Ser43) (PE) Antibody - LS-C281934 - LSBio

<b>CatalogID:</b>	LS-C281934
<b>Target:</b>	v-raf-1 murine leukemia viral oncogene homolog 1 (RAF1)
<b>Synonyms:</b>	RAF1 Antibody, C-Raf Antibody, C-raf-1 Antibody, CRAF Antibody, NS5 Antibody, Raf-1 Antibody, Proto-oncogene c-RAF Antibody, Oncogene RAF1 Antibody, RAF Antibody, Raf1 kinase Antibody
<b>Family / Subfamily:</b>	Protein Kinase / Raf
<b>Host</b>	RAF1 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Phycoerythrin (PE)
<b>Immunogen Species:</b>	RAF1 / RAF antibody was raised against Human
<b>Antigen Type:</b>	Synthetic phospho-peptide
<b>Immunogen:</b>	RAF1 / RAF antibody was raised against kLH-conjugated synthetic phosphopeptide from amino acids around Ser43 of human RAF1.
<b>Specificity:</b>	Human RAF1 / RAF
<b>Epitope:</b>	Ser43
<b>Reactivity:</b>	Human
<b>Purification:</b>	Affinity purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	Western blot, ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences