

RAF1 / RAF Rabbit anti-Human Polyclonal (Ser338) (PE) Antibody - LS-C281998 - LSBio	
CatalogID:	LS-C281998
Target:	v-raf-1 murine leukemia viral oncogene homolog 1 (RAF1)
Synonyms:	RAF1 Antibody, C-Raf Antibody, C-raf-1 Antibody, CRAF Antibody, NS5 Antibody, Raf-1 Antibody, Proto-oncogene c-RAF Antibody, Oncogene RAF1 Antibody, RAF Antibody, Raf1 kinase Antibody
Family / Subfamily:	Protein Kinase / Raf
Host	RAF1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	RAF1 / RAF antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	RAF1 / RAF antibody was raised against kLH-conjugated synthetic phosphopeptide from amino acids around Ser338 of human RAF1.
Specificity:	Human RAF1 / RAF
Epitope:	Ser338
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/30/2014	
© 2014 LifeSpan BioSciences	