

Recoverin Rabbit anti-Human Polyclonal (Internal) (FITC) Antibody - LS-C282091 - LSBio	
CatalogID:	LS-C282091
Target:	recoverin (RCVRN)
Synonyms:	RCVRN Antibody, Recoverin Antibody, Protein CAR Antibody, RCV1 Antibody
Host	RCVRN antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Recoverin antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	Recoverin antibody was raised against kLH-conjugated synthetic peptide from internal region of human Recoverin.
Specificity:	Human Recoverin
Epitope:	Internal
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only Not for resale without prior written consent from LifeSpan BioSciences, Inc. Created on 9/30/2014 © 2014 LifeSpan BioSciences	