

RIPK3 / RIP3 Rabbit anti-Human Polyclonal (C-Terminus) (FITC) Antibody - LS-C282309 - LSBio	
<b>CatalogID:</b>	LS-C282309
<b>Target:</b>	receptor-interacting serine-threonine kinase 3 (RIPK3)
<b>Synonyms:</b>	RIPK3 Antibody, Receptor interacting protein 3 Antibody, RIP-like protein kinase 3 Antibody, RIP3 Antibody, RIP-3 Antibody, Receptor-interacting protein 3 Antibody
<b>Family / Subfamily:</b>	Protein Kinase / RIP
<b>Host</b>	RIPK3 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	RIPK3 / RIP3 antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	RIPK3 / RIP3 antibody was raised against kLH-conjugated synthetic peptide from the C-terminus of human RIPK3.
<b>Specificity:</b>	Human RIPK3 / RIP3
<b>Epitope:</b>	C-Terminus
<b>Reactivity:</b>	Human, Mouse
<b>Purification:</b>	Affinity purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/30/2014	
© 2014 LifeSpan BioSciences	