

C-Kit / CD117 Rabbit anti-Human Polyclonal (C-Terminus) (PE) Antibody - LS-C282542 - LSBio	
CatalogID:	LS-C282542
Target:	v-kit Hardy-Zuckerman 4 feline sarcoma viral oncogene homolog
Synonyms:	KIT Antibody, C-Kit Antibody, CD117 Antibody, PBT Antibody, Piebald trait Antibody, Piebald trait protein Antibody, Proto-oncogene c-Kit Antibody, SCFR Antibody, Soluble KIT variant 1 Antibody, Tyrosine-protein kinase Kit Antibody, p145 c-kit Antibody, CD117 antigen Antibody
Family / Subfamily:	Protein Kinase / PDGF Receptor
Host	KIT antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	C-Kit / CD117 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	C-Kit / CD117 antibody was raised against kLH-conjugated synthetic peptide from the C-terminus of human KIT.
Specificity:	Human C-Kit / CD117
Epitope:	C-Terminus
Reactivity:	Human, Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ICC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/30/2014	
© 2014 LifeSpan BioSciences	