

SUMO Rabbit anti-Drosophila Polyclonal (C-Terminus) (PE) Antibody - LS-C282596 - LSBio	
CatalogID:	LS-C282596
Target:	Small Ubiquitin-Like Modifier (SUMO)
Host	Small Ubiquitin-Like Modifier (SUMO) antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	SUMO antibody was raised against Drosophila
Antigen Type:	Synthetic peptide
Immunogen:	SUMO antibody was raised against kLH-conjugated synthetic peptide from the C-terminus of Drosophila SUMO.
Specificity:	Drosophila SUMO.
Epitope:	C-Terminus
Reactivity:	Drosophila, Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/30/2014	
© 2014 LifeSpan BioSciences	