

SRC Rabbit anti-Human Polyclonal (Tyr419) (PE) Antibody - LS-C283024 - LSBio

<b>CatalogID:</b>	LS-C283024
<b>Target:</b>	v-src avian sarcoma (Schmidt-Ruppin A-2) viral oncogene homolog
<b>Synonyms:</b>	SRC Antibody, ASV Antibody, Pp60c-src Antibody, Tyrosine kinase pp60c-src Antibody, Tyrosine-protein kinase SRC-1 Antibody, v-src Antibody, C-SRC Antibody, p60-Src Antibody, Proto-oncogene c-Src Antibody, SRC1 Antibody
<b>Family / Subfamily:</b>	Protein Kinase / Src
<b>Host</b>	SRC antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Phycoerythrin (PE)
<b>Immunogen Species:</b>	SRC antibody was raised against Human
<b>Antigen Type:</b>	Synthetic phospho-peptide
<b>Immunogen:</b>	SRC antibody was raised against kLH-conjugated synthetic phosphopeptide from amino acids around Tyr419 of human SRC.
<b>Specificity:</b>	Human SRC
<b>Epitope:</b>	Tyr419
<b>Reactivity:</b>	Human
<b>Purification:</b>	Affinity purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences