

SRC Rabbit anti-Human Polyclonal (N-Terminus) (FITC) Antibody - LS-C283028 - LSBio

CatalogID:	LS-C283028
Target:	v-src avian sarcoma (Schmidt-Ruppin A-2) viral oncogene homolog
Synonyms:	SRC Antibody, ASV Antibody, Pp60c-src Antibody, Tyrosine kinase pp60c-src Antibody, Tyrosine-protein kinase SRC-1 Antibody, v-src Antibody, C-SRC Antibody, p60-Src Antibody, Proto-oncogene c-Src Antibody, SRC1 Antibody
Family / Subfamily:	Protein Kinase / Src
Host	SRC antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	SRC antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	SRC antibody was raised against kLH-conjugated synthetic peptide from the N-terminal region of human SRC.
Specificity:	Human SRC
Epitope:	N-Terminus
Reactivity:	Human, Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences