

STAT3 Rabbit anti-Human Polyclonal (Tyr705) (FITC) Antibody - LS-C283092 - LSBio

CatalogID:	LS-C283092
Target:	signal transducer and activator of transcription 3 (acute-phase response factor) (STAT3)
Synonyms:	STAT3 Antibody, Acute-phase response factor Antibody, APRF Antibody, DNA-binding protein APRF Antibody, HIES Antibody
Host	STAT3 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	STAT3 antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	STAT3 antibody was raised against kLH-conjugated synthetic phosphopeptide from amino acids around Tyr705 of human STAT3.
Specificity:	Human STAT3.
Epitope:	Tyr705
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences