

STAT5A Rabbit anti-Human Polyclonal (Ser726) (FITC) Antibody - LS-C283096 - LSBio

CatalogID:	LS-C283096
Target:	signal transducer and activator of transcription 5A (STAT5A)
Synonyms:	STAT5A Antibody, MGF Antibody, STAT5 Antibody
Host	STAT5A antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	STAT5A antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	STAT5A antibody was raised against kLH-conjugated synthetic phosphopeptide from amino acids around Ser726 of human STAT5a.
Specificity:	Human STAT5A / STAT5.
Epitope:	Ser726
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences