

SYK Rabbit anti-Human Polyclonal (Tyr525/526) (PE) Antibody - LS-C283284 - LSBio

CatalogID:	LS-C283284
Target:	spleen tyrosine kinase (SYK)
Synonyms:	SYK Antibody, p72syk Antibody, Spleen tyrosine kinase Antibody, Syk protein-tyrosine kinase Antibody, Tyrosine-protein kinase SYK Antibody, p72-Syk Antibody
Family / Subfamily:	Protein Kinase / Syk/ZAP70
Host	SYK antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	SYK antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	SYK antibody was raised against kLH-conjugated synthetic phosphopeptide from amino acids around Tyr525/526 of human Syk.
Specificity:	Human Syk.
Epitope:	Tyr525/526
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences