

TAZ Rabbit anti-Human Polyclonal (N-Terminus) (FITC) Antibody - LS-C283390 - LSBio

CatalogID:	LS-C283390
Target:	tafazzin (TAZ)
Synonyms:	TAZ Antibody, BTHS Antibody, EFE Antibody, Endocardial fibroelastosis 2 Antibody, EFE2 Antibody, LVNCX Antibody, Protein G4.5 Antibody, Tafazzin Antibody, CMD3A Antibody, G4.5 Antibody, TAZ1 Antibody, XAP-2 Antibody
Host	TAZ antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	TAZ antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	TAZ antibody was raised against kLH-conjugated synthetic peptide from the N-terminal region of human TAZ.
Specificity:	Human TAZ
Epitope:	N-Terminus
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences