

MAPT / TAU Rabbit anti-Human Polyclonal (C-Terminus) (FITC) Antibody - LS-C283391 - LSBio	
CatalogID:	LS-C283391
Target:	microtubule-associated protein tau (MAPT)
Synonyms:	MAPT Antibody, FTDP-17 Antibody, MSTD Antibody, Neurofibrillary tangle protein Antibody, PHF-tau Antibody, MTBT2 Antibody, PPND Antibody, DDPAC Antibody, MAPTL Antibody, MTBT1 Antibody, Paired helical filament-tau Antibody, TAU Antibody
Host	MAPT antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	MAPT / TAU antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	MAPT / TAU antibody was raised against kLH-conjugated synthetic peptide from the C-terminus of human TAU.
Specificity:	Human MAPT / Tau
Epitope:	C-Terminus
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/30/2014 © 2014 LifeSpan BioSciences	