

TERT / Telomerase Rabbit anti-Human Polyclonal (Ser1125) (FITC) Antibody - LS-C283465 - LSBio

CatalogID:	LS-C283465
Target:	telomerase reverse transcriptase (TERT)
Synonyms:	TERT Antibody, DKCB4 Antibody, EST2 Antibody, HEST2 Antibody, PFBMFT1 Antibody, Telomerase Antibody, TP2 Antibody, DKCA2 Antibody, HTRT Antibody, TCS1 Antibody, Telomerase catalytic subunit Antibody, TRT Antibody
Host	TERT antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	TERT / Telomerase antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	TERT / Telomerase antibody was raised against kLH-conjugated synthetic phosphopeptide from amino acids around Ser1125 of human TERT.
Specificity:	Human TERT / Telomerase
Epitope:	Ser1125
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences