

TLR1 Rabbit anti-Mouse Polyclonal (C-Terminus) (FITC) Antibody - LS-C283650 - LSBio

CatalogID:	LS-C283650
Target:	toll-like receptor 1 (TLR1)
Synonyms:	TLR1 Antibody, CD281 antigen Antibody, KIAA0012 Antibody, Rsc786 Antibody, Toll-like receptor 1 Antibody, TIL. LPRS5 Antibody, TIL Antibody, CD281 Antibody
Family / Subfamily:	Toll-like Receptor / not assigned-Toll-like Receptor
Host	TLR1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	TLR1 antibody was raised against Mouse
Antigen Type:	Synthetic peptide
Immunogen:	TLR1 antibody was raised against kLH-conjugated synthetic peptide from the C-terminus of mouse TLR1.
Specificity:	Mouse TLR1.
Epitope:	C-Terminus
Reactivity:	Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences